

NEW!

LM

LMArte™ Solo

feel the difference

INSTRUMENT DUO FOR SIMPLIFYING YOUR DAILY CASES

Streamline your restorative instrumentation with this duo that is all you need for your daily anterior and posterior direct composite restoration cases.

WHY?

These instruments are ideal for simplifying your daily cases: with just two instruments, you can perform most of your high-end anterior and posterior restorations.

WHAT?

LM-Arte™ Solo instruments simplify everything from modeling to ergonomics and even hygiene, with fewer instruments needed on the tray.



Sharp tip
Precision modeling
Micromodeling
Limit marking
Scratching

Thin blade
Proximal modeling

Convex profile
One-point modeling
Create depressions
Create buccal profile
Create cervical profile
Access difficult areas

Flat profile
Big excess removal
Medium excess removal

Flat end
Plugger

Straight edges
Scraping

Flat profile
Flattening
Uniform

Angled edge
Definition
Scraping

Angulation
Access for deep cavities

Semi-rounded edges
Slope forming
Excess removal
Margin removal



4-5 mm bulkfill safety area
Probe for bulk fill materials

Tapered profile
Fine excess removal
Slope modeling
Secondary anatomy

Sharpened edges
Cut, split, separate
Marginal ridges

Rounded tip
Delicate delimitation
Cervical contouring

Sharp tip
Precision modeling
Micromodeling
Limit marking
Scratching
Flow distribution
Bubble elimination



ONE INSTRUMENT FOR DIRECT VENEERING

LM-Arte™ Solo Anterior LM 431-433 ES

Instrument for high-end esthetic anterior restorations. Rounded and straight spatulas are designed for modeling anterior teeth, especially labial surfaces.

ONE INSTRUMENT FOR POSTERIORES

LM-Arte™ Solo Posterior LM 434-438 ES

Instrument for high-end posterior composite restorations. Plugger end is designed for bulk filling and composite material modeling. The sharp end enables contouring the occlusal surface anatomy.



Learn more at LM™ website www.lm-dental.com/LM-Arte