

HOW COULD YOU BENEFIT FROM LM Dental Tracking System™?



Are you doing fact based budgeting and purchasing?



ASSET MANAGEMENT

- Process optimization
- Easy consumption recording
- Daily clinic management
- Budgeting & long term planning



COST REDUCTION

- Decrease losses
- Improve material flow
- Optimize stock
- Avoid time consuming manual controls



How are you improving the threshold for patient safety?



PATIENT SAFETY

- Infection control documentation
- Automated checks and alerts
- Material traceability
- Records of maintenance and sterilization procedures

Are you able to document the instruments use of your students?



EDUCATION SUPPORT

- Validate competences in simulation exercises
- Feedback to students
- Documented exercises



Check your potential profit at
www.dentaltracking.com/profit-calculator

Scan the QR
to see the
DTSTM™ animation



ARE YOU USING DATA

to benchmark operation performance and patient outcomes for constant improvement of best practices?

RFID identification tags

The LM-ErgoSense hand instruments and other dental materials with identification tags can be scanned with DTS™ readers for reliable tracking through their entire life cycle - from storage and maintenance to clinical and patient care.



DTS™ readers

The DTS™ readers (either multi-read or single-read) can be placed in any environment. They notify users of errors, such as the use of unsterilized instruments or expired material.

DTS™ server software

The DTS™ server is used to record, track and verify the cycles of activities done to tagged materials.

Analysis

The software provides easy analytical reports about the items and activity cycles registered in the system.

Action

The information in the reports can improve and ease asset management and documentation, and most importantly increase cost efficiency and patient safety.

Interested to know more?

Contact: dts.sales@lm-dental.com | www.dentaltracking.com

LM Dental

LM-Instruments Oy | Norrbyn rantatie 8, FI-21600 Parainen, Finland
Tel. +358 2 4546 400 | info@lm-dental.com | www.lm-dental.com