

# LM-PerioKit

LM

feel the  
difference



STYLE ITALIANO™  
MADE IN ITALY

# LM-PerioKit

## LM™ Periodontal Kit powered by STYLEITALIANO™

LM-PerioKit by Styleitaliano™ is a set of LM™ periodontal instruments chosen by Dr. Giano Ricci.

LM™ periodontal instruments are ultimate in metallurgy. The high wear resistance and hardness of the blade metal makes them efficient, long-lasting and comfortable. The ergonomic design with a silicone handle gives the instruments a non-slip and natural grip. Good tactile sensitivity and easy handling of the instruments guarantee the best clinical results.



LM-PerioKit LM 6842 ES



**North Carolina Probe**  
LM 54B ES

For detecting calculus, cavities and abnormalities on the tooth surfaces as well as probing the depth of pockets. With 5-10-15 mm markings.



**Furcation Probe**  
LM 20B-21B ES

For examination of furcation lesions from different angles in both the upper and lower jaws.



**McCall 17S-18S**  
LM 279-280 ES

For the removal of sub- and supra-gingival calculus. Excellent in tight interdental spaces.



**Concavity Diamond File SV 5-6**  
LM 255-256D ES

For acute furcation lesions. Diamond coated round curved toe.

## “The LM-PerioKit by StyleItaliano™ contains the essentials for the best possible therapy”

Instruments have been selected in the spirit of StyleItaliano™ whose philosophy is to make treatment feasible, teachable, repeatable. Periodontal therapy too can follow this approach. It's my feeling that, in the best interest of his patients, every practitioner should be able to properly execute non surgical therapy, crown lengthening procedures, free gingival grafts, simpler surgical root coverages.”

**Giano Ricci**  
**MD. DDS. MScD.**



## LMErgoSense® handle

“Feels great in the hand”

Found to have the highest rated ergonomics in studies by the Finnish Institute of Occupational Health: design prevents unnecessary clenching of the hand and makes working pleasant.



**Mini Gracey 7/8**  
LM 207-208M ES

For premolars and molars. For buccal, mesial and lingual surfaces. Great for deep and narrow pockets.



**Mini Gracey 11/12**  
LM 211-212M ES

For premolars' and molars' mesial surfaces. Great for deep and narrow pockets.



**Mini Gracey 13/14**  
LM 213-214M ES

For premolars' and molars' distal surfaces. Great for deep and narrow pockets.



**Implant Mini Universal curette**  
LM 283-284MTI EM

Titanium tips for sensitive but still effective implant maintenance. Universal design works well in tight interdental spaces. With replaceable LM-ErgoMix™ tips.

# LM Dental™



LM-Instruments Oy

Norrbyn rantatie 8 | FI-21600 Parainen | Finland

Tel. +358 2 4546 400

info@lm-dental.com | [www.lm-dental.com](http://www.lm-dental.com)



## LM-PerioKit by Styleitaliano™

Styleitaliano™ is a study group of passionate clinicians and researchers specialized in aesthetic restorative dentistry and periodontal treatments.

More information at [www.styleitaliano.org](http://www.styleitaliano.org)  
and [www.facebook.com/styleitaliano](https://www.facebook.com/styleitaliano).

