

EXCHANGEABLE TIPS

LM-ErgoMix

LM-ErgoMix hand instrument series offers the dental team a specially economical, ecological and ergonomic alternative. Whenever there is a need for a tip replacement the lock grip is opened, a new tip is inserted and the lock grip is closed again. No extra tools are needed.

Use LM-ErgoMix instruments in the normal way, clean and sterilize it like any other instrument. The sealed construction guarantees high level hygiene.

Benefits of retipping:

- Always sharp instruments
- Minimizes waste
- Enables personalized tip combinations.



TITANIUM TIPS FOR IMPLANTS

A great instrument series for implant maintenance. The ErgoMix structure ensures excellent tactile sensitivity and optimal rigidity.



Tips are made of sensitive but still effective titanium, which will not harm the abutments.

Implant Mini Universal curette

LMErgoMix LM 283-284MTI EM

- Delicate working end
- Especially for tight interdental spaces



Implant Mini Gracey 1/2

LMErgoMix LM 201-202MTI EM

- For incisors and premolars



Implant Mini Gracey 11/12

LMErgoMix LM 211-212MTI EM

- For premolars' and molars' mesial surfaces



Implant Mini Gracey 13/14

LMErgoMix LM 213-214MTI EM

- For premolars' and molars' distal surfaces

LM-ERGOMIX SPARE PARTS

Spare handle

LM-ErgoMix handle dark grey LM 9011 EM

Spare lock grips

Orange LM 9020 EM

Light blue LM 9021 EM

Grey LM 9022 EM

Yellow LM 9023 EM

Dark grey LM 9024 EM

Red LM 9025 EM

Spare titanium tips

Implant Mini Universal Curettel L LM 283MTI EM

Implant Mini Universal Curette R LM 284MTI EM

Implant Mini Gracey 1 LM 201MTI EM

Implant Mini Gracey 2 LM 202MTI EM

Implant Mini Gracey 11 LM 211MTI EM

Implant Mini Gracey 12 LM 212MTI EM

Implant Mini Gracey 13 LM 213MTI EM

Implant Mini Gracey 14 LM 214MTI EM