

LM DarkDiamond™

LM

feel the
difference

PLEASANT AND PRECISE
SCULPTING OF COMPOSITE



LM

feel the
difference

Diamond-like feel & precision

Restorative instruments with LM Dark Diamond™ non-stick coating are ideal for placing and contouring of composite resins. The diamond-like coating reduces stickiness and allows pleasant and precise restorative instrumentation.



LM

feel the
difference

- Exclusive diamond-like coating for non-stick placement and sculpting of composite resin
- Not sensitive to scratches due to hard surface
- Smooth surface for easy cleaning
- Dark-colored coating maximizes contrast against fillings and minimizes reflections
- Ergonomic LM-ErgoSense® technology delivers excellent precision and feel

[Click here for a video about LM Dark Diamond™](#)



LM Dark Diamond™

Instrument range

Burnisher

- LM Dark Diamond™ non-stick coating
- For packing and modeling composite

LM ErgoSense®

LM 31-32DD ES



Burnisher-Plugger, ball-ended 1.5-2.5 mm

- LM Dark Diamond™ non-stick coating
- For packing and modeling composite

LM ErgoSense®

LM 471-473DD ES



LMDarkDiamond™

Instrument range

Modeller MH

- LM Dark Diamond™ non-stick coating
- For packing large proximal fillings
- Cylindrical working end will not get stuck in composite during light-curing

LMErgoSense®

LM 708-709DD ES




Modeller MH Mini

- LM Dark Diamond™ non-stick coating
- For packing small and medium proximal fillings
- Cylindrical working end will not get stuck in composite during light-curing

LMErgoSense®

LM 707-708DD ES



LM Dark Diamond™


Instrument range

Packer-Modeller IO

- LM Dark Diamond™ non-stick coating
- For packing and modeling composites
- Conical working end optimal for shaping occlusal surfaces

LM ErgoSense®

LM 482-483DD ES



Plastic instrument, anterior

- LM Dark Diamond™ non-stick coating
- For placing, packing and contouring composite

LM ErgoSense®

LM 43-47DD ES



LMfeel the
difference

LM Dark Diamond™

Instrument range

Plastic instrument, posterior

- LM Dark Diamond™ non-stick coating
- For placing, packing and contouring composite

LM ErgoSense®

LM 48-702DD ES



Mini Spatula

- LM Dark Diamond™ non-stick coating
- For placing and contouring composite

LM ErgoSense®

LM 444-445DD ES



LM Dark Diamond™ instruments can be cleaned and sterilized like any other LM™ instrument.

LM Dental™



LM-Instruments Oy

Norrbyn rantatie 8 | FI-21600 Parainen | Finland

Tel. +358 2 4546 400

info@lm-dental.com | www.lm-dental.com

